


JAIMY

Motorized Articulated Instrument

Open your **M**ind, Extend your **H**ands

**Unlimited
rotation**



5 mm

Ergonomic

Articulated



Dexterity


5 mm

Autoclavable



jaimy
the only **5mm** motorized
articulated instrument

The first motorized articulating laparoscopic instrument with **iD-intelligent Dexterity**, designed to enable surgeons to overcome the unique challenges presented by single incision and conventional laparoscopic surgery.



***iD-Ring Control:**
control movements from your fingertip with the fist intelligent one finger directed instrument. Intuitive and instinctive, iD-Ring Control puts all instrument movements at the tip of your finger!

- ✓ **Full articulation** that maps the surgeon's hand motion
- ✓ **Reusable** for improved economics and cost containment
- ✓ Motorized distal movements for enhanced **dexterity** and **ergonomy**
 - Motorization provides fine movements and instrument rigidity for demanding applications
 - **Total flexion** for unparalleled access to reach over, under, and around structures
 - **Unlimited 360° rotation** of the end effector at any bending angle for **easy suturing** and tissue manipulation
- ✓ Instrument movement and control functions are easily and ergonomically accessed with **one hand** (fits left or right hand)
- ✓ Provides **intelligent dexterity** with seven degrees of freedom for improved access, superior tissue manipulation, and easy **triangulation** in even the most challenging cases
- ✓ Benefits a wide variety of **surgical procedures** in different clinical specialties, including urology, bariatric surgery, general surgery, gynecology and cardiothoracic surgery



Non contractual document